


BIO 201 Lab 1

Experiments 1, 2, 3

Professor Diane Hilker



Overview

- I. Exp. 1: Introduction to the Microscope
 - II. Exp. 2: Survey of Microbes
 - III. Exp. 3: Collection of Microbes
- 

I. Exp. 1: Intro. to the Microscope

► **Purpose:** To review the use & care of the compound light microscope

Ocular lens (eyepiece)
Remagnifies the image formed by the objective lens

Body tube Transmits the image from the objective lens to the ocular lens

Arm

Objective lenses
Primary lenses that magnify the specimen

Stage Holds the microscope slide in position

Condenser Focuses light through specimen

Diaphragm Controls the amount of light entering the condenser

Illuminator Light source

Coarse focusing knob

Base

Fine focusing knob

(a) Principal parts and functions

Copyright © 2010 Pearson Education, Inc.




Fig. 3.1 Textbook

I. Exp. 1: Intro. to the Microscope

- ▶ **Compound Binocular Light Microscope**
 - **(2) Sources of Magnification:**
 - Eyepiece or Ocular (10x)
 - Objectives (4):
 - Scanning Power: 4 x
 - Low Power: 10x
 - High Power: 40x
 - Oil Immersion: 100x
- ▶ **Parfocal:** ability to go from one objective to another with minimal focusing

I. Exp. 1: Intro. to the Microscope

▶ Total Magnification: TM

- $TM = \text{Magnification of Eyepiece} \times \text{Magnification of Objective}$

	Eyepiece Magnification	Objective Magnification	TM
Scanning	10X	4X	40X
Low	10X	10X	100X
High Dry	10X	40X	400X
Oil Immersion	10X	100X	1000X

I. Exp. 1: Intro. to the Microscope

► Resolution or Resolving Power (RP)

- Ability to distinguish detail clearly
- To be able to tell 2 points as separate points and not one point


- $RP = \frac{\text{Wavelength of Light}}{2 \times \text{Numerical Aperture}} = \frac{\lambda}{2NA}$

I. Exp. 1: Intro. to the Microscope

- ▶ Wavelength of light (nm)

Red light = 700 nm

Blue light = 400 nm



- ▶ $RP_{\text{red}} = \frac{700 \text{ nm}}{2(1)} = 350 \text{ nm}$

- ▶ $RP_{\text{blue}} = \frac{400 \text{ nm}}{2(1)} = 200 \text{ nm}$

- ▶ Lower the resolution, better the clarity
- ▶ Blue filter provides the best resolution with a halogen light bulb
- ▶ Blue filter NOT needed with microscopes that have a LED light bulb

I. Exp. 1: Intro. to the Microscope

- ▶ **Numerical Aperture (NA):** describes the cone of light that enters the lens so as to see fine detail. Two things make up NA:
 - **Angular Aperture:** angle of light as it goes through the lenses & filters of the condenser & into the objective (Constant)
 - **Refractive Index:** how light travels through a medium
 - Refractive Index of Air = 1.0
 - Refractive Index of Oil = 1.5

I. Exp. 1: Intro. to the Microscope

▶ $RP_{\text{air}} = \frac{400 \text{ nm}}{2(1.0)} = 200 \text{ nm}$

▶ $RP_{\text{oil}} = \frac{400 \text{ nm}}{2(1.5)} = 133 \text{ nm}$

- ▶ **Better resolution
with oil**




Overview

- I. Exp. 1: Introduction to the Microscope
- II. Exp. 2: Survey of Microbes
- III. Exp. 3: Collection of Microbes

II. Exp. 2: Survey of Microbes


- ▶ **Purpose:** To become familiar with using a microscope & to view various microbes
- **Wet Mount:** observing living cells
 - Focus on edge of coverslip
 - Scanning—dim light using diaphragm
 - Move toward center of slide
 - Observe under Low & High Powers
 - Slides will dry out quickly



II. Exp. 2: Survey of Microbes

► 4 Slides: Largest to smallest microorganisms

1. Pond Water: algae–much variation



II. Exp. 2: Survey of Microbes

2. **Protozoa:** single celled eukaryotic microbes that move by different methods that belong to the Protista kingdom.

- **Pseudopods:** false feet

Amoeba



- **Cilia**

Paramecium



- **Flagella**


Euglena



II. Exp. 2: Survey of Microbes

3. Yeast: single celled eukaryotic microbes that belong to the Fungi kingdom


- Ovoid & irregular
- Budding: method of reproduction
- Brownian movement
- Smaller than protozoa
- Larger than bacteria



II. Exp. 2: Survey of Microbes

4. Bacteria (Hay infusion): single celled prokaryotic microbes that belong to the Monera kingdom.

- Must view under 400x
- Very small
- Motile & non-motile
- Looks like specks of sand
- Hard to discern shape
- Smaller than yeast
& protozoa
- Protozoa may be present in the sample



Overview

- I. Exp. 1: Introduction to the Microscope
- II. Exp. 2: Survey of Microbes
- III. Exp. 3: Collection of Microbes

III. Exp. 3: Collection of Microbes

- ▶ **Purpose:** To collect and grow microbes from the environment for observation
 - Procedure to be described by lab instructor