

BIO 201 Lab 2

Experiment 3

Professor Diane Hilker

Overview

I. **Exp. 3: Collection of Microbes**

1. Observe different types of Microbial Colonies
2. Identification of Molds
3. Isolation of Bacteria

Overview

I. **Exp. 3: Collection of Microbes**

1. Observe different types of Microbial Colonies

2. Identification of Molds

3. Isolation of Bacteria


I. Exp. 3: Collection of Microbes

1. Microbial Colonies


- **Colony:** a visible mass of microbial cells originating from one cell.
- **(3) Types**
 - Large, fuzzy, hairy, 3D, growing upward & touching the lid, various colors–**MOLD**
 - Small, creamy, moist, circular, various colors–**BACTERIA**
 - Medium, dry, crusty, white–beige–**YEAST**

I. Exp. 3: Collection of Microbes

1. Microbial Colonies



Mold Colonies



Bacterial Colonies



Yeast Colonies

I. Exp. 3: Collection of Microbes

▶ Culture Media Used

◦ Potato Dextrose Agar (PDA)

- Supports mainly mold growth
- pH 5.2–acidic
- High in carbohydrates

◦ Nutrient Agar (NA)

- Supports mainly bacterial growth
- pH 7.0–neutral
- High in proteins

Overview

I. **Exp. 3: Collection of Microbes**

1. Observe different types of Microbial Colonies

2. Identification of Molds

4. Isolation of Bacteria


I. Exp. 3: Collection of Microbes

▶ Molds

Vegetative Structures: obtain nutrients


- Absorb nutrients through cell wall
- Can't identify a mold based on vegetative structure

I. Exp. 3: Collection of Microbes

- **Thallus:** body of mold consisting of filaments
 - **Hyphae or hypha:** filaments–multicellular
 - Can be very long; elongate at the tips
 - **Septa or septum:** cross–walls
 - **Coenocytic hyphae:** no cross–walls
 - **Mycelium:** filamentous mass visible to the eye
- 

I. Exp. 3: Collection of Microbes

Fig. 12.2 Textbook



(a) Septate hypha

(b) Coenocytic hypha

(c) Growth of a hypha from a spore

I. Exp. 3: Collection of Microbes

- ▶ **Molds–Reproductive Structures: Spores**
 - Molds identified based on **type of spores**
 - 2 Types of Spores:
 - **Sexual Spores:** genetic exchange between 2 parents (meiosis)
 - Not as common in nature
 - To be discussed in lecture
 - **Asexual Spores:** no genetic exchange (mitosis)
 - More common in nature
 - To be discussed in lab

I. Exp. 3: Collection of Microbes

▶ Asexual Spores: 2 Types




I. Conidiospores or conidia: 2 types based on size of spore

1. Microconidia

- Conidiophore: supporting structure
- Holds conidia
- Examples: *Penicillium sp.* and *Aspergillus sp.*



I. Exp. 3: Collection of Microbes

▶ *Penicillium sp.*



I. Exp. 3: Collection of Microbes

▶ *Aspergillus sp.*



I. Exp. 3: Collection of Microbes

▶ Asexual Spores: 2 Types

I. Conidiospores or conidia: 2 types based on size of spore

2. Macroconidia: much larger than microconidia

- Examples: *Alternaria*, *Stemphyllium*, *Stachybotrys*, *Curvularia*, *Fusarium*, and *Microsporium*

I. Exp. 3: Collection of Microbes

▶ Macroconidia

Alternaria



Curvulvaria



Stemphyllium



Fusarium




Microsporium



I. Exp. 3: Collection of Microbes

- ▶ Macroconidia
Stachybotrys



I. Exp. 3: Collection of Microbes



▶ Asexual Spores: 2 Types

II. Sporangiospores

- Sporangium: sac
- Sporangiophore: supporting structure
 - Holds sporangiospores
- Examples: *Rhizopus sp.* and *Mucor sp.*

I. Exp. 3: Collection of Microbes

▶ *Rhizopus sp.* and *Mucor sp.*



I. Exp. 3: Collection of Microbes

- ▶ How to make a slide under a Biological Safety Cabinet (BSC)
 - Choose a sporulating mold colony
 - Place 2 drops of ethanol on slide
 - Aseptically remove a small but visible piece of the mold using acceptable tools
 - Add 1 drop of **Lactophenol Cotton Blue**
 - Cover with cover slip
 - Observe under Scan (dim light), Low and High Power if needed

I. Exp. 3: Collection of Microbes

▶ Observing Prepared Slides

- Observe prepared mold slides using Scan, Low and High Powers
 - *Rhizopus sp.* – sporangiospores
 - *Penicillium sp.* – conidiospores (microconidia)
 - *Aspergillus sp.* – conidiospores (microconidia)
- Observe hyphae
- Can you tell the difference between sporangiospores and conidiospores?

Overview

I. **Exp. 3: Collection of Microbes**

1. Observe different types of Microbial Colonies


2. Identification of Molds

3. Isolation of Bacteria

I. Exp. 3: Collection of Microbes

▶ Isolation of Bacteria

- **Pure Culture:** 1 type of microbe; to get alone
- Procedure called **Streaking** or **Streaking for Isolation**
- **T-Streak Method:** to separate individual colonies



- To be demonstrated by instructor