


BIO 201 Lab 9
Experiments 13,14, 17, & 18
Results

Professor Diane Hilker



Overview

- I. **Exp. 13: Evaluation Antimicrobial Agents**
 - II. **Exp. 14: Evaluation of Antibiotics**
 - III. **Exp. 17: Skin Flora–*Staphylococcus***
 - IV. **Exp. 18: Throat Flora**
- 


I. Exp. 13: Antimicrobial Agents

- ▶ Fill in Table 7 with Ring Sizes (0 to a 5)




- ▶ Which disinfectant worked the best?
- ▶ Which disinfectant worked the least?

Overview

- I. Exp. 13: Evaluation Antimicrobial Agents
 - II. Exp. 14: Evaluation of Antibiotics
 - III. Exp. 17: Skin Flora–*Staphylococcus*
 - IV. Exp. 18: Throat Flora
- 

II. Exp. 14: Antibiotics

- ▶ Fill in Table 8 with Ring sizes (0 to 5)



Copyright © 2007 Pearson Education, Inc., publishing as Benjamin Cummings.

- ▶ Which is a broad spectrum antibiotic?
- ▶ Which is a narrow spectrum antibiotic?

Overview

- I. Exp. 13: Evaluation Antimicrobial Agents
- II. Exp. 14: Evaluation of Antibiotics
- III. Exp. 17: Skin Flora–*Staphylococcus*
- IV. Exp. 18: Throat Flora

III. Exp. 17: Skin Flora

3 Types of Culture Media Used

1. Enrichment Media : **Nutrient Agar (NA)**
 2. Selective Media–selecting for halophiles
 - **Both *S. aureus* & *epidermidis* will grow
- **Staph 110**
 - **Mannitol Salt Agar (MSA)**

III. Exp. 17: Skin Flora

3 Types of Culture Media Used

3. Differential Media


- **Mannitol Salt Agar (MSA)**

- *S. aureus*: plate yellow

Why? Ferments mannitol,

pH drops, phenol red turns yellow.

- *S. epidermidis*: plate stays pink because it doesn't ferment mannitol



Overview

- I. Exp. 13: Evaluation Antimicrobial Agents
- II. Exp. 14: Evaluation of Antibiotics
- III. Exp. 17: Skin Flora–*Staphylococcus*
- IV. Exp. 18: Throat Flora

IV. Exp. 18: Throat Flora


- ▶ **Purpose:** To isolate & examine microbes obtained from the throat & to observe the three different types of hemolytic reactions.
- ▶ **Differential Culture Media**
 - **Blood Agar Plates (BAP):** 3 different types of hemolysis

IV. Exp. 18: Throat Flora

Blood Agar Plates: 3 Types of Hemolytic Reactions

1. Gamma Hemolysis: γ

- No breakdown of rbc around the colony
- Rbc's intact



2. Alpha Hemolysis: α


- Partial breakdown of rbc around the colony
- Rbc's are yellowish, greenish, brownish



IV. Exp. 18: Throat Flora


3. Beta Hemolysis: β

- Complete breakdown of rbc around the colony
- Rbc's are broken down and clear
- Indicates a “higher degree of pathogenicity”



IV. Exp. 18: Throat Flora

- ▶ Strep. Throat
 - Beta-hemolytic *Streptococcus pyogenes*



III. Exp 17

- ▶ Save either a Staph 110 plate or MSA plate that shows medium sized beige colonies
 - Assume the colony is *Staphylococcus*
 - Do the following:
 1. **Gram Stain:** Gram pos. cocci in clusters
 - Heat fix & Gram Stain per lab manual
 2. **Catalase Test** (Refer to Exp. 15 in the Lab Manual)

III. Exp 17

▶ Catalase Test:




◦ Procedure:

1. Using aseptic technique add a small amt. of the colony to the center of a slide
2. **Do NOT add water or heat fix**
3. Add 1–2 drops of hydrogen peroxide
4. Look immediately for bubbling or fizzing as a result of the O_2 gas being given off

III. Exp 17

▶ Catalase Test:

- *Staphylococcus*: Catalase positive–bubbles
- *Streptococcus*: Catalase negative–no bubbles



BIO 201 Lab 9

Experiments 19 & 20

Professor Diane Hilker




Overview

- I. **Exp. 19: Water Analysis for Fecal Contamination**
- II. **Exp. 20: Quantitative Analysis of H₂O**

I. Exp. 19: Water Analysis for Feces

- ▶ **Purpose:** To learn & perform the 3 stage standard H₂O analysis test for fecal contamination
 - **PRLB Tube:** Phenol Red Lactose Broth with a Durham Tube (collects gas)
 - Follow Instructor's directions.



II. Exp. 20: Quantitative Analysis of Water

- ▶ **Purpose:** To test the hypothesis that potable H₂O may still contain bacteria and be safe to drink. To compare the number of bacteria in tap vs. well H₂O.
 - ▶ Tap Water vs. Well Water
 - ▶ Follow Instructor's Directions
- 